

ROOSEVELT SPECIFICATIONS

Minimum weight: 35 lbs.
Maximum weight: 115 lbs.

Minimum height: 33.5 inches
Maximum height: 62 inches

Shoulder width maximum: 15.5 inches
Thigh width maximum: 14.5 inches

Shoulder strap heights:

- 12.5"
- 15"
- 17.25"
- 19.75"
- 22.25"

Padded seat extenders provide optimum seat depth for each leg. Any combination of depths is available.

| Left Leg | Right Leg |
|----------|-----------|
| 12.5" | 12.5" |
| 15.5" | 15.5" |
| 17.0" | 17.0" |

Height from seat bottom to seat top is 30 inches.

Our "EZ-Tether" system makes tethering unnecessary in most cases, simplifying the installation in older vehicles.

Ask us how we can customize your Roosevelt!

As always, we accept returns and charge NO restocking fees!

See DEMO VIDEOS:
www.youtube.com/user/MerrittManufacturing



Designed and manufactured in the United States by Merritt Manufacturing, Inc., The Roosevelt™ conforms to all applicable Federal motor vehicle safety standards.

©2016 Merritt Manufacturing, Inc. All Rights Reserved

MERRITT MANUFACTURING, INC.

PO Box 17152
Indianapolis, IN 46217
Phone: 317-409-0148
Fax: 317-893-2567
Email: info@merrittcarseat.com

www.MerrittCarSeat.com
Check us out on social media!



THE ROOSEVELT™

A car seat specially designed for children with special health care needs.



OPTIONAL FEATURES

OTHER OPTIONS:
Padded Seat Extenders, Extra Long Harness

"EZ-UP" HEADREST SYSTEM



Developed for children with poor head control, the EZ-Up system enables the child to maintain a natural head position. Use with the EZ-Up Cap™ to restrict forward motion.

"QUICK-CHANGE" INCONTINENT COVER



Eliminate time consuming removal of cover, harness, etc. Just remove the leak proof incontinent cover, wash and replace in no time.

SCOLIOSIS HARNESS



Independent bilateral hip adjusters and extra foam pads allow therapists to optimally position the child by filling any gaps behind the child's back.

"EZ-GUARD" BUCKLE COVER



Designed to prevent the occupant from releasing the car seat harness buckle.

CHEST CLIP GUARD



The chest clip guard prevents the occupant from pushing down on or unfastening the chest clip.

"STAY-PUT" POMMEL



The pommel (abductor) stays firmly in place while in use. Pivots away for easy transfer.