

P.O. Number:

Date:

## 2022 Churchill Order Form

The Churchill-Pediatric Positioning Vehicle Restraint System and Accessories

## See Important Statement on Last Page About Installing the Churchill\*\*\*\*

The Churchill is available in three fabric color options. To select a color, simply fill in the space at the end of the part number in the table with the color option abbreviation from the following list.

Sapphire Blue – S Total Camo – TC Beige Naugahyde - N

## **Billing Address**

Contact Name:	
Organization:	
Street Address:	
City, State:	
ZIP Code:	
Phone:	Fax:
E-Mail:	

## **Shipping Address (If Different)**

Contact Name:	
Organization:	
Street Address:	
City, State:	
ZIP Code:	
Phone:	Fax:
E-Mail:	

Qty	Part Number	Description	Unit Price	Line Total
	2000CS15	The Churchill Standard* with Vest and 15 Inch Seat Depth – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Vest is included in price. Please choose vest size below.		
	2000CS18	The Churchill Standard* with Vest and 18 Inch Seat Depth – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Vest is included in price. Please choose vest size below.		
	2000CH15	The Churchill with Positioning Harness*** and 15 Inch Seat Depth (FOR USE WITH OUR ESCAPE PROOF PRODUCTS) – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Harness is included in price. NOTE: Add the escape proof Buckle Guard and Chest Clip Guard listed below for the proficient escape artist.		

1			I	
	2000CH18	The Churchill with Positioning Harness*** and 18 Inch Seat Depth (FOR USE WITH OUR ESCAPE PROOF PRODUCTS) – a pediatric positioning vehicle restraint system (Weight range: 65 to 175 lbs.; Height Range: 48 to 72 inches). One Positioning Harness is included in price. NOTE: Add the escape proof Buckle Guard and Chest Clip Guard listed below for the proficient escape artist.		
	2000CVS	Small Positioning Vest - for chest circumference between 25 and 32 inches and belly circumference less than 32 inches. <b>**</b> Price charged only when ordering vest separately.		
	2000CVM	Medium Positioning Vest - for chest circumference between 30 and 37 inches and belly circumference less than 37 inches. <b>**</b> Price charged only when ordering vest separately.		
	2000CVL	Large Positioning Vest - for chest circumference between 35 and 44 inches and belly circumference less than 44 inches. <b>**</b> Price charged only when ordering vest separately.		
	2000VB	Add Side Release Buckles to one of the vests above. (Side release buckles can help stop an occupant from escaping the vest. For the proficient escape artist, rather than the vest with buckles, we suggest selecting the harness version with Buckle Guard and Chest Clip Guard listed below).		
		Optional Accessories		
	2000BG	Buckle Guard – ONLY FOR USE WITH 2000CH15-S or 2000CH18-S		
	2000CCG	Chest Clip Guard – ONLY FOR USE WITH 2000CH15-S or 2000CH18-S		
	2000EUHR	EZ-Up Head Rest – Supports head to limit side to side motion. To limit forward motion, choose a cap size from below: The cap is designed to Velcro to the head rest.		
	See note to the right regarding the three cap options.	Note: There are three options for cap color. Sapphire Blue (S), Total Camo (TC), or Khaki (K). The Khaki option is the coolest option for hot temperatures because it is single layer. Simply replace the asterisk in the part numbers below with "S" for Sapphire Blue, "TC" for Total Camo, or "K" for Khaki.		
	1000EUCS	EZ-Up Cap — Size Small - (Head Circumference: 17.5 — 19.5 inches)		
	1000EUCM	EZ-Up Cap — Size Medium - (Head Circumference: 19.5 — 21.5 inches)		
	1000EUCL	EZ-Up Cap — Size Large - (Head Circumference: 21 — 23.25 inches)		
	2000FA	Foam Abductor		
	2000HICP	Hip Inflection Comfort Pommel		
	2000HIW	Hip Inflection Wedge		

2000lC-15	Quick-Change Incontinent Cover for 15 Inch Seat Depth Base		
2000IC-18	Quick-Change Incontinent Cover for 18 Inch Seat Depth Base		
2000ADD	Adductor Strap		
2000AB	Attachment Bracket Kit (for Vehicles without LATCH) – <b>See</b> Important Note Below****		
-		Sub Total**	
	Indiana Residents Only; Not Applicable to other States	Sales Tax (7%)	
		TOTAL	

* The "Churchill Standard" includes the following	***The "Churchill with Harness" includes the	
features:	following features:	
<ul> <li>Base</li> <li>Lower Anchor Attachment Straps</li> <li>Positioning Harness for Vest Option</li> <li>Anti-Slouching Leg Strap Assembly</li> <li>**Positioning Vest (Small, Medium, or Large)</li> <li>Foam Padding</li> <li>Fabric Cover</li> <li>Tether Extension Assembly</li> <li>Belt Positioning Guides</li> </ul> ** One Vest is included in price of "Standard Churchill with Vest." Retail Price of additional vests is \$275.00 each.	<ul> <li>Base</li> <li>Lower Anchor Attachment Straps</li> <li>Positioning Harness for 5-Point Harness Option</li> <li>5-Point Harness</li> <li>Foam Padding</li> <li>Fabric Cover</li> <li>Tether Extension Assembly</li> <li>Belt Positioning Guides</li> </ul>	

\*\*\*\* The Attachment Bracket Kit is designed for use in vehicles that do not have LATCH. The Churchill must be installed with LATCH or with the Attachment Bracket kit. All Vehicles 2003 and newer are equipped with LATCH in at least two rear seating positions, so there should be no installation issues for these vehicles and the Attachment Bracket will not be necessary. Many 2000 - 2002 models also have LATCH (See your vehicle owner's manual to verify if LATCH is present).

For vehicles without LATCH, the attachment bracket can <u>only</u> be used on vehicle seats that you can slide a ruler, or similar object, from the back side of the vehicle seat through the crack in the vehicle seat (where the seat bottom and seat back come together) through to the front side of the seat. This can usually be done in mini-vans and SUVs. Generally, the Attachment Bracket cannot be installed in the back seat of sedans. Installation of the Attachment Bracket Kit does require that a hole be drilled through the vehicle's floor board.

Please call Customer Service at 317-409-0148 with any questions regarding the installation of the Churchill.

Please place orders by one of the following methods:
E-Mail: order@merrittcarseat.com

	Phone: 317-409-0148 Fax: 317-893-2567
	Thanks for Your Business!
** Shipping charges – The cost of shipping "Best Way" will be added to Purchase Orders.	Merritt Manufacturing, Inc., 400 N. Forsythe Street, Unit 1 Franklin, IN 46131
	www.merrittcarseat.com